

Certificate of Accreditation

Blood Bank Computer Systems, Inc. (BBCS) certifies
that the following is true and correct:

BBCS and Terumo BCT (Vendor)

worked together to ensure:

The Terumo BCT Vista Information System, Version 1.1 developed by the vendor can interface with the ABO Wheels Application, Version 1.1.0 when the Vista Interface developed by BBCS is used.

BBCS tested over 40 transactions, during which data that was sent on the Vista Interface was verified by BBCS and was able to be properly retrieved, transmitted, and loaded for display using the Vista Information System, Interface, and the ABO Wheels Application together.

After passing the test, BBCS certifies that the Terumo BCT Vista Information System is compatible with the ABO Wheels Application, Version 1.1.0. When used together pursuant to the respective operating instructions for each, the Vista System and the ABO Wheels Application are capable of retrieving, transmitting and loading for display the intended and correct data.

This certification is limited to the software versions identified herein, and is not applicable to any modified or later version tested, except as expressly provided herein.

BBCS Accreditation Program: Device Certification



Certificate of Accreditation

Outbound Data Elements

Outbound data elements are sent to the Vista System in preparation to receive data from the collection device. Because it is not end result data, outbound data elements are not subject to the same level of testing as inbound data elements.

Data Elements

| |
|-------------------|
| Donation Site |
| Appearance Date |
| Appearance Time |
| Donor Number |
| Donor SSN |
| Donor Last Name |
| Donor First Name |
| Donor Middle Name |
| Donor Race |
| Sex |
| Blood Type |
| DOB |
| Address |
| City |
| State |
| Zip Code |
| Day Phone |
| Night Phone |

Data Elements

| |
|---|
| Email |
| Donation Type |
| Donation Product |
| Pre-procedure Platelet Count |
| CMV Status* |
| Hematocrit [◊] |
| Hemoglobin [◊] |
| CuSO4 Flag |
| Donor Height |
| Donor Temperature |
| Pulse Rate |
| Systolic |
| Diastolic |
| Platelet Product Eligibility Calculation |
| Plasma Product Eligibility Calculation |
| Red Cell Product Eligibility Calculation |
| Whole Blood Product Eligibility Calculation |
| Comments |

* It is a known product issue that the Terumo BCT Vista Information System, Version 1.1 will translate blank/null values in the CMV field to NEG, creating potential for a false-negative to be displayed within Vista. Please note that the CMV field in ABO Wheels is not impacted by this issue at any point.

[◊] Only one test type (HGB or HCT) may be used.



Certificate of Accreditation

Certified Inbound Data Elements

| Data Elements | Included in Interface |
|-----------------------------|-----------------------|
| Pre-platelet Count | Certified |
| Product Code | Certified |
| Arm Inspection (Either Arm) | Certified |
| TU Code | Certified |
| Draw Start Date/Time | Certified |

| Data Elements | Included in Interface |
|-------------------------|-----------------------|
| Draw Stop Date/Time | Certified |
| Phlebotomist ID – Start | Certified |
| Phlebotomist ID – Stop | Certified |
| Plasma Loss | Certified |
| Red Blood Cell Loss | Certified |

